G-aenial LIGHT-CURED RADIOPAQUE COMPOSITE RESTORATIVE For use by a dental professional in the recommended indications.

- RECOMMENDED INDICATIONS A. G-aenial ANTERIOR 1. Direct restorative for Class III, IV, V cavities. 2. Direct restorative for wedge-shaped defects and root surface cavities. 3. Direct restorative for veneers and diastema closure. B. G-aenial POSTERIOR 1. Direct restorative for Class I and II cavities.

- CONTRAINDICATIONS 1. Fuly capping. 2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician. DIRECTIONS FOR USE 1. Shade selection. Clean the tooth with pumice and water. Shade selection should be made prior to isolation. Select the appropriate G-aenial shades by referring to the G-aenial shade guide.

- 2. Cavity preparation. Prepare cavity using standard techniques. Dry by gentle blowing with oil free air. Note: For full capping, use calcium hydroxide. 3. Bonding Treatment For bonding G-aenial to enamel and / or dentin, use a light-cured bonding system such as GC G-aenial Bond, GC G-BOND™, GC Fuji BOND LC or GC UniFil® Bond (Fig. 1). Follow manufacturer's instructions. 4. Placement of G-aenial 1) Dispensing from a Syringe Insert the G-aenial Unitip into the Unitip APPLIER or equivalent. Remove the cap and extrude material directly into the prepared cavity. 2) Dispensing from a syringe Remove syringe cap and dispense material onto a mixing pad. Place the material into the cavity using a suitable placement instrument. After dispensing, screw syringe plunger anticlockwise by a half to turn to release residual pressure inside the syringe. Replace cap immediately after use.

- Note: 1. Basically, material can be applied in a single layer to achieve aesthetic restorations using Standard shades. For details, refer to the Clinical Hints. 2. Material may be hard to extrude immediately after removing from cold storage. Prior to use, leave to stand for a few minutes at room temperature. 3. After dispensing, avoid too long exposure to ambient light. Ambient light can shorten the manipulation time.

- Clinical Hints a. In the case of small cavities Restoring a small cavity using a single shade layering technique. In most cases the use of one Standard shade alone will be sufficient. In cases where a higher degree of translucency is needed, one of the Outside special shades can be selected. See also Examples of Clinical Applications. b. In the case of large cavities In most cases a multi shade layering technique will give the best aesthetic results. To block out shine throughs from the oral cavity or shade discoloration dentin, select an appropriate inside special shade and continue to build up with a Standard shade. To make a restoration more life like e.g. to copy age related changes in appearance, the final layer should consist of an Outside special shade. See also Examples of Clinical Applications and / or consult the Shade Composition Chart.

- 2. Posterior cavities a. In the case of small cavities Restoring a small cavity using a single shade layering technique. In most cases the use of one Standard shade alone will be sufficient. In cases where a higher translucency is needed, one of the Outside special shades can be selected. See also Examples of Clinical Applications. b. In the case of large cavities In most cases a multi shade layering technique will give the best aesthetic results. To block out shine throughs from the oral cavity or shade discoloration dentin, select an appropriate inside special shade and continue to build up with a Standard shade. To make a restoration more life like e.g. to copy age related changes in appearance, the final layer should consist of an Outside special shade. See also Examples of Clinical Applications and / or consult the Shade Composition Chart.

- 3. Contouring before Light Curing Contour using standard techniques. Light cure G-aenial using a light curing unit (Fig. 3). Keep light guide as close as possible to the surface. Refer to the following chart for Irradiation Time and Effective Depth of Cure.

G-aenial ANTERIOR: Irradiation Time and Effective Depth of Cure

Table with columns for Irradiation time (3 sec, 6 sec, 10 sec, 20 sec, 40 sec) and Effective Depth of Cure (3.0 mm, 3.5 mm, 2.5 mm, 2.0 mm, 1.5 mm).

G-aenial POSTERIOR: Irradiation Time and Effective Depth of Cure

Table with columns for Irradiation time (3 sec, 6 sec, 10 sec, 20 sec, 40 sec) and Effective Depth of Cure (3.0 mm, 3.5 mm, 2.5 mm, 2.0 mm, 1.5 mm).

STORAGE

Store in a cool and dark place (4 – 25°C / 39.2 – 77.0°F) away from high temperatures or direct sunlight.

PACKAGES

- 1. Quick Start Kit - 7 syringes (1 each in 7 anterior shades: A1, A2, A3, B2, AE, IE, JE) (2.7mL per syringe) - Shade guide

- 2. REFILL - 1 syringe (in 28 shades) (2.7mL per syringe) Note: Weight per syringe: 4.7g for anterior syringe, 4.0g for TE shade, 5.5g for posterior syringe

- III. Option 1. Unitip APPLIER 2. G-aenial SHADE GUIDE 3. Mixing pad (No. 14B)

CAUTION

- 1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. 2. In case of contact with eyes, flush immediately with water and seek medical attention. 3. Take care to avoid ingestion of the material. 4. Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent possible sensitivity. 5. For infection control reasons, Unitips are for single use only. 6. Wear protective eye glasses during light curing. 7. When polishing the polymerized material, use dust collector and wear a dust mask to avoid inhalation of cutting dust. 8. Do not mix with other similar products. 9. Avoid getting material on clothing. 10. In case of contact with unhardened areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing. 11. Do not use G-aenial in combination with eugenol containing materials as eugenol may hinder G-aenial from setting. 12. All shades except for the outside special shade TE are radioopaque.

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