**everX Posterior™**  
FIBER REINFORCED COMPOSITE FOR DENTIN REPLACEMENT

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

- Cavities including 3 surfaces or more
- Cavities with missing cusps
- Deep cavities (including class I, II and endodontically treated teeth)
- Cavities after amalgam replacement
- Cavities where onlays & inlays would also be indicated

Note:  
everX Posterior should always be covered with a layer of light-cured universal restorative composite, for sufficient wear resistance

**CONTRAINDICATIONS**

- Do not use for pulp capping.
- At least one horizontal dimension of the cavity should exceed 3 mm.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- The product is not suitable for patients with a history of hypersensitivity to methacrylate monomers.
- Do not use as final approximal and surface composite layer.

**DIRECTIONS FOR USE**

- Shade Selection

Clean the tooth with pumice and water. Shade selection of the final composite enamel layer should be made prior to isolation.

- Cavity Preparation (Fig.1)

Prepare cavity using standard techniques. Dry by gently blowing with oil free air. The use of rubber dam is recommended to isolate prepared tooth from contamination with saliva, blood or sulcus fluid.

Note:  
For pulp capping, use calcium hydroxide.

- Bonding treatment (Fig.2)

For bonding everX Posterior to enamel and/or dentine, use a light-cure/cure bonding system such as GC G-aenial Bond, GC G-BOND. Follow manufacturer's instructions.

- Placement of everX Posterior (Fig.3) – see note for Class II and large cavities

Insert everX Posterior into the Unitip APPLIER or equivalent. Remove the cap, and extrude everX Posterior material directly into the cavity.

Pack the everX Posterior in the cavity, using suitable instruments, until it has been adapted tightly against the cavity floor and walls. To prevent premature curing, avoid bright light in the work area. When placing everX Posterior material in the cavity, ensure that there is enough space for the overlaying composite on all surfaces of the restoration (refer to point 6).

- Class II and large cavities (Fig.3a, 3b)

a. Use a suitable matrix system to establish the proximal contour and contact area.

b. When placing everX Posterior in Class II or larger cavity, build first the missing walls by using light-curing restorative composite before placing everX Posterior. The wall should be thick enough to withstand the application pressure of everX Posterior.

Note:  
Material may be hard to extrude immediately after removing from cold storage. Prior to use, leave to stand for a few minutes at normal room temperature.

**Clinical Hint**

Optionally, apply first a thin layer of flowable composite (e.g., G-aenial Flo, Gradia Direct Flo) to the cavity floor before the application of everX Posterior. Place the everX Posterior on top of the flowable composite and pack it in the cavity.

- Light Curing (Fig.4)

Light cure everX Posterior using a light curing unit. Keep light guide as close as possible to the surface. Refer to the following chart for Irradiation Time and Effective Depth of Cure.

Curing device	Irradiation Time	Depth of cure
Plasma arc (2000 mW/cm²)	9 sec.	4 mm
High Power LED Light (>1200 mW/cm²)	10 sec.	4 mm
Halogène/ Normal LED Light (700 mW/cm²)	20 sec.	4 mm

Note:  
Lower light intensity may cause insufficient curing or discoloration of the material.

- Placement of the overlaying composite (Fig.5)

Fill the remaining cavity space with a wear-resistant and polishable restorative composite, such as GC G-aenial Posterior or GC Kalore. The composite layer should be of 1-2 mm on the occlusal surface.

Light-cure as instructed by the manufacturer.

- Finishing and Polishing

Check and adjust the occlusion. Finish and polish using standard techniques.

**SHADE**

everX Posterior is available in one universal shade

**STORAGE**

Store in a cool and dark place (4 - 25°C / 39.2 - 77.0°F) away from high temperatures or direct sunlight.

(Shelf life : 2 years from date of manufacture)

**PACKAGES**

15 unitips (0,25g / 0,13mL)

**CAUTION**

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Take care to avoid ingestion of the material.
- Do not mix with other similar products.
- Do not mix with other similar products.
- For infection control reasons, Unitips are for single use only.
- Wear rubber-free plastic or rubber gloves during operation
- direct contact with air inhibited resin layers in order to prevent possible sensitivity.
- Wear protective eye glasses during light curing.
- When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
- Avoid getting material on clothing.
- In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing.
- Do not use everX Posterior in combination with eugenol containing material which may hinder everX Posterior from setting.
- Do not wipe the unitip with excessive alcohol as this may cause contamination of the paste.
- Dispose of all wastes according to local regulations.

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