Prior to use, carefully read the instructions for use.

OSTRON 100

ACRYLIC TRAY AND BASEPLATE RESIN

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

Construction of individual impression trays and baseplates.

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE

A. INDIVIDUAL TRAY

- PREPARATION
 - Take an impression for making a study model.
 - To create space for the impression material and for covering undercuts, apply a suitable spacer to the model.
- MIXING
 - Standard powder/liquid ratio is 10 g / 4.5 g and is obtained by one scoop powder (13 g) and a full measure of liquid (5.8 g). Two scoops of powder are sufficient for an upper arch, and one and a half scoops for lower.
 - Mix them in the bowl for 30 seconds.
- FORMING
 - Roll out the mixture uniformly to a thickness of 1.5 to 2 mm.
 - Place the flattened mixture over the spacer, form a tray and trim off excess.
 - Make a handle with residual mixture before hardening. This can be done by applying liquid to both tray and handle.
 - When retention holes are needed, the use of amalgam-carrier before hardening or round bur after hardening is recommended.
- TRIMMING
 - After Ostron becomes hard, remove it from the model and brush the inside clean. If a wax spacer is used, it must be washed out with hot water containing detergent.
 - Smooth the edges using a suitable instrument.
- SETTING
 - Apply the impression material adhesive to the inside and edge of the tray, and leave until it dries.
 - Load the impression material into the tray and set it into the mouth.

B. BASE PLATE

- Soak the model in water or apply a separator (ACROSEP) to the mucosal surface of the model.
- After mixing, roll out the mixture and form the base plate on the model.
- Trim off the excess on the model using a wax carver before hardening.

STORAGE

Recommended for optimal performance, store in a cool place away from direct sunlight.

PACKAGES

- Powder: 500 g tin, 1 kg tin, 5 kg tin, 10 kg tin (Europe: 1 kg only)
- Liquid: 250 g (270 mL) bottle, 2.5 kg (2.7 L) tin (Europe: 250 g only)

SHADES

Transparent blue, Blue, Transparent pink (Europe: Transparent blue, Transparent pink only).

CAUTION

- The liquid is volatile and highly flammable. Close the container immediately after use and keep away from sources of ignition.
- Avoid direct contact between skin and uncured resin. It is recommended to wear rubber or plastic gloves.
- Avoid contact with skin or eyes. In case of contact with eyes, flush immediately with water and seek medical attention.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com>
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

Undesired effects - Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

Last revised: 01/2020

**'GC'**

MANUFACTURED by
GC DENTAL PRODUCTS CORP.
2-285 Tomimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED by
GC CORPORATION
76-1 Hiasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Luven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza, Singapore 529541
TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200

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