G-CEM™ CAPSULE

SELF ADHESIVE LUTING CEMENT IN CAPSULES

For use only by a dental professional in the recommended indications.

- RECOMMENDED INDICATIONS**
- Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges (max. 2 abutments)
 - Cementation of metal-ceramic, fiber posts, and cast post and cores, and Sitaflauten.

- CONTRAINDICATIONS**
- Pulp capping.
 - Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE

Powder / Liquido Ratio (g/g)	0,35 / 0,15
Mixing Time (sec.)	10"
Working Time (23°C/73°F) (min., sec.)	2'30"

A. CEMENTATION TECHNIQUE FOR CERAMIC, INDIRECT COMPOSITES and metal based restorations

- 1. TOOTH PREPARATION**
 - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
 - Clean the prepared tooth with pumice and water.
 - Rinse thoroughly with water. Dry by gently blowing with oil free air.

- 2. RESTORATION PREPARATION**
 - Prepare the restorations pre-treated and handled according to the manufacturer's instructions.

- 3. CAPSULE ACTIVATION AND MIXING**
 - Before use, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
 - Place the capsule into the plunger until it is fully depressed with the main body and hold it down for 2 seconds (Fig. 2).

Mixing time chart

Mixer Speed	Mixing Time
Standard mixer (approx. 4.000rpm)	10 seconds
High-speed mixer (over 4.500rpm)	10 seconds

- 4. DISPENSING**
 - Immediately remove the mixed capsule from the mixer and attach it to the GC CAPSULE APPLIER.
 - Make two clicks to the capsule (Fig. 3) to adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.

- 5. CEMENTATION**
 - Coat the internal surface of the restoration with sufficient cement and seal immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73°F). Higher temperatures will shorten working time. Retrioperator of G-CEM CAPSULE is recommended to extend the working time.

- 6. EXCESS CEMENT REMOVAL**
 - To remove excess cement, light cure using a visible light curing unit (Halogen LED) for 2 minutes. Do not allow the cement to self-cure until the excess cement feels rubbery.
 - Remove excess cement by scraping (Fig. 7).
 - Light cure each surface for 10 seconds (single surface from occlusal: 10 seconds, any other surfaces: 10 seconds each) or allow the cement to self cure for 4 minutes after setting of the restoration (Fig. 9).

- 7. FINISH AND POLISH**
 - Finish and polish margins with appropriate instruments.

B. CEMENTATION TECHNIQUE FOR metal, fiberglass and ceramic posts

- 1. TOOTH PREPARATION**
 - Prepare the post space according to manufacturer's instructions. The post space should be cleaned, rinsed and dried thoroughly using paper points.
 - The post space should be dry.

- 2. POST PREPARATION**
 - Make sure that the post is pretreated and handled according to the manufacturer's instructions.

- 3. CAPSULE ACTIVATION AND MIXING**
 - Refer to A.3.

- 4. DISPENSING & CEMENTATION**
 - Refer to A.4.

- 5. EXCESS CEMENT REMOVAL**
 - Refer to A.6.

- 6. STORAGE**
 - Recommended for optimal performance, store the capsules in the original aluminum foil in a cool and dark place (4-25°C) (39,2-77°F).

- SHADES**
 - A2 (Vita® shade), A03 (Opak), B01 (Opak) and Translucunt

- PACKAGE**
 1. Single shade package: 50 capsules
 2. Assorted package: 50 capsules (20 of each A2 and Translucunt, 5 each of A03 and B01)

- AVERAGE CONTENTS per capsule:** 0,35 g powder and 0,15 g (0,14 mL) liquid
- Minimum net volume of mixed cement per capsule:** 0,17 mL
- GC CAPSULE APPLIER:** 1 piece
- GC Capsule Elongation Tip:** 50 pieces

- CAUTION**
 - In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.
 - In case of contact with eyes, flush immediately with water and seek medical attention.
 - If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks.
 - Do not use G-CEM CAPSULE in combination with eugenol containing materials as eugenol may hinder G-CEM CAPSULE from setting or bonding properly.
 - Use a protective light shield or similar protective eye wear during light curing.
 - This product is not indicated for filling or core build-up.
 - When using GC Capsule Elongation Tip with a cotton soaked with alcohol, prevent the alcohol from entering the capsule nozzle as this will affect the setting of the resin.
 - In rare cases the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
 - Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.


Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for the Americas: <http://www.gcamerica.com>

They can also be obtained from your supplier.

Undesired effects: Warning!

If you become aware of any kind of undesired effect, reaction or similar events experienced by you or this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/rectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental in this way you will contribute to improve the safety of this product.

Last revised: 06/2019

MANUFACTURED BY
GC CORPORATION
7-1, Higashi-Azabu 1-chome, Minato-ku, Tokyo 106-8585, Japan
RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Anso, B. 80003 U.S.A.
TEL: +1 303 426 6000

GC CORPORATION
7-1 Higashi-Azabu 1-chome, Minato-ku, Tokyo 106-8585, Japan
RESPONSIBLE MANUFACTURER IN EUROPE:
GC AMERICA INC.
3737 West 127th Street, Anso, B. 80003 U.S.A.
TEL: +1 303 426 6000

GC SOUTH AMERICA
Rua Helena 9, 506, Tupy, São Paulo, SP, BRASIL
CNPJ: 06.278.959/0001-01
CNPJ: 06.278.959/0001-01
CNPJ: 06.278.959/0001-01

GC ASIA DENTAL PTE. LTD.
11 Tampines Street 1, #01-01 Tampines Plaza, Singapore 520541
TEL: +65 6546 7388

GC AUSTRALASIA DENTAL PTY. LTD.
111 Springvale Road, Springvale, VIC 3171, Australia
TEL: +61 3 9301 8000

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- 6. ZEMENTAZIONE**
 - Coat the internal surface of the restoration with sufficient cement and seal immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73°F). Higher temperatures will shorten working time. Retrioperator of G-CEM CAPSULE is recommended to extend the working time.

- 7. EXCESS CEMENT REMOVAL**
 - To remove excess cement, light cure using a visible light curing unit (Halogen LED) for 2 minutes. Do not allow the cement to self-cure until the excess cement feels rubbery.
 - Remove excess cement by scraping (Fig. 7).
 - Light cure each surface for 10 seconds (single surface from occlusal: 10 seconds, any other surfaces: 10 seconds each) or allow the cement to self cure for 4 minutes after setting of the restoration (Fig. 9).

- 8. FINISH AND POLISH**
 - Finish and polish margins with appropriate instruments.

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