Prior to use, carefully read the instructions for use. EN

**GC Fuji TEMP LT™**

GLASS IONOMER PROVISIONAL LUTING CEMENT

For use only by a dental professional in the indications for use.

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Paste B: Contrast agent, water, silicon dioxide, polyacrylic acid, glycerine, polyacrylic carboxylic acid, stabilizer, paraben.

**IDENTIFICATION OF PARTS (Fig. 1)**

- A. Paste Pak Cartridge  
1. Cover 2. Inner plunger 3. Cartridge tip  
B. Paste Pak Dispenser  
4. Piston 5. Adjusting guide 6. Lever 7. Sliding block 8. Piston release

**PASTE PAK CARTRIDGE LOADING**

To retract the piston fully, pull the piston release and at the same time pull the sliding block (on the underside of the dispenser) with fingers.

Load the cartridge into the Paste Pak Dispenser (Fig. 3). Ensure the mark (▲) on the cartridge lines up with the mark (▲) on the Paste Pak Dispenser as illustrated. Then turn the cartridge fully to the right position.

After loading, push the piston release forward until it stops (Fig. 4).

- a) Push the piston release forward.  
b) Do not push the upper side.

NOTE:  
To remove the cartridge, retract the piston fully and turn the cartridge back to the marked positions.**DIRECTIONS FOR USE**

Paste / Pasto Ratio (g/g)	0.7 / 1.0
Mixing Time (sec.)	10'-15'
Working Time (23°C/73°F) (min., sec.)	2'00"
Finishing Time (37°C/99°F, after removing excess) (min., sec.)	3'00"

1. **ABUTMENT PREPARATION**

1) Clean the tooth/abutment surfaces with pumice and water. For pulp capping, use calcium hydroxide.

2) Rinse abutment with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESCALICATE. Prepared surfaces may cause difficulty in cleaning the residual cement after debonding. NOTE:

a. When a moderate retention is desired and especially when the fit of restoration is tight, a thin layer of separation medium like petroleum jelly can be applied on the abutment surface. Prior to final, permanent cementation, the petroleum jelly can be removed by the use of petroleum jelly hand setting or bonding of the permanent cement.

b. The use of desensitizer may lower the bonding of GC Fuji TEMP LT.

It is recommended to increase the surface of the restoration when increased retention is desired. For sandblasting procedure, refer to the respective manufacturer's instructions. Before cementation, apply a thin layer of petroleum jelly to the outer marginal surface of the restoration. This will ease the clean-up of excess cement.

**3. DISPENSING**

1) Turn over the cartridge cover.

2) Depress the lever to dispense the required amounts of paste onto the mixing pad.

a) Before dispensing, flatten the material from the cartridge onto the mixing pad.

b) The adjusting guide can be moved forwards or backwards to change the quantity of material to be dispensed (Fig. 5).

1. Position of adjusting guide for one way.

2. Position of adjusting guide for two ways.

c) The paste may be dispensed in the shape of a rod, but this does not affect the quality of the paste.

3) Depress the lever while moving the cartridge and dispenser backward (Fig. 6), level off and cut the extruded material from the cartridge tip by using the cartridge and dispenser to an upright position on the mixing pad (Fig. 7).

4) When the lever is released after dispensing, residual pastes at the cartridge tip are drawn into the cartridge. Clean the cartridge tip and replace the cartridge cover.

NOTE:

a. The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the underside of the dispenser.

4. **MIXING**

First spread the yellow paste along the two pastes together with a spatula. Spread patiently on the mixing pad and spatulate vigorously and thoroughly with lapsing strokes for 10 seconds. Be careful not to incorporate air bubbles.

NOTE:

a. Add a larger amount of paste requires a longer mixing time (15 seconds).

b. Do not mix GC Fuji TEMP LT with other materials.

5. **CEMENTATION**

1) Apply a thin layer of mixed GC Fuji TEMP LT on the internal rim of the restoration. The working time is 2 minutes from start of mixing at 23°C (73°F). Higher temperature will shorten working time.

2) Seal the restoration. Maintain moderate pressure.

a) Before sealing, allow the cement to self-cure until excess cement feels rubbery (approx. 1 minute).

b) After sealing, allow the restoration to set.

c) When the lever is released after dispensing, residual pastes at the cartridge tip are drawn into the cartridge. Clean the cartridge tip and replace the cartridge cover.

NOTE:

a. The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the underside of the dispenser.

4. **RESTORATION REMOVAL**

1) Use a suitable crown lifter and dental forceps to remove previously cemented restorations. Avoid any force that could damage the preparation or restoration.

2) After removal of the restoration, residual cement remaining on the preparation and/or internal surfaces of the restoration, may be removed with a ultrasonic scalpel, explorer or other suitable hand instrument such as a sharp scalpel. Gently remove the cement from the preparation and the preparation.

3) Treat the preparation and proceed to final cementation according to the cement manufacturer's instructions or repeat provisional cementation with GC Fuji TEMP LT Cement by following initial cementation procedure noted at the beginning of this instruction for use document.

**CLEANING / STERILIZATION OF Paste Pak Dispenser**

For cleaning and sterilization of Paste Pak Dispenser, refer to the respective instructions for use.

**SHADE**

Universal

Recommended for optimal performance, store at room temperature away from heat, moisture and direct sunlight (4-25°C) (39.2-77°F).

**PACKAGE**

1) GC Fuji TEMP LT Cartridge (13.3 g / 7.2 mL) (2), mixing pad (1)

2) Paste Pak Dispenser (1)

**CAUTION**

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and / or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue. When using the product, please follow the instructions for use.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. Avoid getting material on clothes.

4. In rare cases the product may cause sensitivity in some people. If any such reactions occur, discontinue the use of the product and refer to a physician.

5. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

http://www.gceurope.com

or for the Americas:

http://www.gcamerica.com

They can also be obtained from your supplier.

**CLEANING AND DISINFECTION**

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection immediately after use inspect device and label for use if the device is damaged.

DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Use a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUADAMED database (<https://ec.europa.eu/tools/euadamed>) or contact us at Regulatory.gce@gc.dental

Undesired Effects - Response:

In case of contact with any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in the instructions for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

http://www.gceurope.com/vigilance/gc.dental

In this way you will contribute to improve the safety of this product.

UK Responsible Person

GC UNITED KINGDOM Ltd.

Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom

Last revised: 12/2023

R Only

CE 2797

UK CA 0086

/GC/

MANUFACTURED by

GC CORPORATION

76-1 Hakenoma-cho, Habata-ku, Tokyo 174-8585, Japan

RESPONSIBLE MANUFACTURER IN CANADA:

377 West 17th Street, Aliso, IL 6003 U.S.A.

DISTRIBUTED BY

GC CORPORATION

76-1 Hakenoma-cho, Habata-ku, Tokyo 174-8585, Japan

EU: GC AMERICA, Inc.

Research Park Haarleweg 1240, Internevelaan 33,

TEL: +45 6546 7588

PRINTED IN JAPAN

Prior to use, carefully read the instructions for use. EN

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**PASTE PAK CARTRIDGE LOADING**

To retract the piston fully, pull the piston release and at the same time pull the sliding block (on the underside of the dispenser) with fingers.

Load the cartridge into the Paste Pak Dispenser (Fig. 3). Ensure the mark (▲) on the cartridge lines up with the mark (▲) on the Paste Pak Dispenser as illustrated. Then turn the cartridge fully to the right position.

After loading, push the piston release forward until it stops (Fig. 4).

- a) Push the piston release forward.  
b) Do not push the upper side.

NOTE:

To remove the cartridge, retract the piston fully and turn the cartridge back to the marked positions.

**DIRECTIONS FOR USE**

Paste / Pasta Ratio (g/g)	0.7 / 1.0
Mixing Time (sec.)	10'-15'
Working time (23°C/73°F) (min., sec.)	2'00"
Finishing Time (37°C/99°F, after removing excess) (min., sec.)	3'00"

1. **ABUTMENT PREPARATION**

1) Clean the tooth/abutment surfaces with pumice and water. For pulp capping, use calcium hydroxide.

2) Rinse abutment with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESCALICATE. Prepared surfaces may cause difficulty in cleaning the residual cement after debonding.

NOTE:

Overmixing may cause difficulty in cleaning the residual cement after debonding.

a. When a moderate retention is desired and especially when the fit of restoration is tight, a thin layer of separation medium like petroleum jelly can be applied on the abutment surface. Prior to final, permanent cementation, the petroleum jelly can be removed by the use of petroleum jelly hand setting or bonding of the permanent cement.

b. The use of desensitizer may lower the bonding of GC Fuji TEMP LT.

It is recommended to increase the surface of the restoration when increased retention is desired. For sandblasting procedure, refer to the respective manufacturer's instructions. Before cementation, apply a thin layer of petroleum jelly to the outer marginal surface of the restoration. This will ease the clean-up of excess cement.

**3. DISPENSING**

1) Turn over the cartridge cover.

2) Depress the lever to dispense the required amounts of paste onto the mixing pad.

a) Before dispensing, flatten the material from the cartridge onto the mixing pad.


b) The adjusting guide can be moved forwards or backwards to change the quantity of material to be dispensed (Fig. 5).

1. Position of adjusting guide for one way.

2. Position of adjusting guide for two ways.

c) The paste may be dispensed in the shape of a rod, but this does not affect the quality of the paste.

3) Depress the lever while moving the cartridge and dispenser backward (Fig. 6), level off and cut the extruded material from the cartridge tip by using the cartridge and dispenser to an upright position on the mixing pad (Fig. 7).

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