
UNIFAST™ III

SELF-CURING ACRYLIC RESIN FOR TEMPORARY INLAYS, CROWNS, BRIDGES AND REPAIRS

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Construction of temporary inlays, crowns and bridges.
- Repair of fractured dentures.

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE

- TEMPORARY INLAYS, CROWNS AND BRIDGES**
 - Direct Method – single shade standard technique
 - Take an impression prior to tooth preparation
 - Coat the prepared teeth and gingival tissue with a thin film of vaseline.
 - Dispense the desired amount of UNIFAST III powder in a rubber cup and add liquid. Quickly mix with a plastic spatula for 10-15 seconds. Immediately pour the mixture into the intended area of the impression. Standard powder to liquid ratio is 1 g of powder to 0.5 mL of liquid (1 scale of powder) to 1 scale of liquid. Powder to liquid ratio is adjustable between 2 g : 1 mL and 2 g : 1.5 mL.
 - Wait until the mixture reaches a soft dough stage.
 - Remove the impression in the mouth with pressure.
- NOTE:** UNIFAST III starts to set 1 minute after mixing. Make sure to finish above steps within this period. If a longer working time is desired, store the powder and liquid in a refrigerator for more than 4 hours.
- Remove the impression from the mouth before the exothermic reaction starts. Approximate setting time is 1 minute and 30 seconds from start of mixing. At this stage, material still has a rubbery consistency.
- Allow the UNIFAST III to bench cure for 10 minutes. The setting time is 3 minutes 10 seconds from start of mixing.
- Remove the temporary appliance from the impression.
- Polish the surface of the temporary appliance with a fine abrasive and establish proper marginal adaptation and occlusion. The margins can be relined by adding a layer of freshly mixed UNIFAST III.

- NOTE:** In case of making a temporary bridge using pre-fabricated temporary crowns as pontics, immerse the bridge in water at 50°C for 5 minutes or leave undisturbed at room temperature for 30 minutes before placing it in the patients' mouth to ensure complete setting.
- After final polishing, lime the inlay/crown/bridge with a temporary cement like FREEGENOL TEMPORARY PAKC.

CLINICAL HINT

The surface of the temporary appliance can be coated with G-COAT PLUS for a high gloss esthetic result. Refer to the respective manufacturer's instructions for use.

- Direct Method – multi shade standard technique
 - Use a tooth matching Vita shade in the technique as described above (Direct method, single shade technique 1-9).
 - Trim the incisal area of the temporary restoration.
 - Apply an incisal enamel or incisal resin to the incisal area, reusing the impression or by applying with a brush (brush-on technique).
 - Continue with Direct method 9).
- Indirect Method – Refer to the respective manufacturer's instructions for use.
 - Fabricate a stone model using a standard technique and take an impression / silicone index.
 - Wax-up restoration using a standard technique and take an impression / silicone index.
 - Coat the stone model with a separator such as ACROSEP or MULTI-SEP and place the fractured denture base on the stone model.
 - Dispense the desired amount of UNIFAST III powder in a rubber cup and add liquid. Quickly mix with a plastic spatula for 10-15 seconds. Immediately pour the mixture into the desired area of the impression/silicone index. Standard powder to liquid ratio is 1 g of powder to 0.5 mL of liquid (1 scale of powder) to 1 scale of liquid. Powder to liquid ratio is adjustable between 2 g : 1 mL and 2 g : 1.5 mL.
 - Remove the model with pressure.
 - Place the impression on the model with pressure.

- NOTE:** UNIFAST III starts to set 1 minute after mixing. Make sure to finish above steps within this period. If a longer working time is desired, store the powder and liquid in a refrigerator for more than 4 hours.
- Remove the temporary appliance from the impression/silicone index.
- Polish the surface of the temporary appliance using a cotton pellet soaked with ethanol.
- Adjust and finish in usual manner. Carefully fit on the model and establish proper marginal adaptation and occlusion. The margins can be relined by adding a layer of freshly mixed UNIFAST III.

- NOTE:** In case of making a temporary bridge using pre-fabricated temporary crowns as pontics, immerse the bridge in water at 50°C for 5 minutes or leave undisturbed at room temperature for 30 minutes before placing it in the patients' mouth to ensure complete setting.

CLINICAL HINT

The surface of the temporary appliance can be coated with OPTIGLAZE for a high gloss esthetic result. Refer to the respective manufacturer's instructions for use.

- After final polishing, lime the inlay/crown/bridge with a temporary cement like FREEGENOL TEMPORARY PAKC.
- Indirect Method – multi shade standard technique
 - Use a tooth matching Vita shade in the technique as described above (Indirect method, 1-8).
 - Trim the incisal area of the temporary restoration.
 - Apply an incisal enamel or incisal resin to the incisal area, reusing the impression/silicone index or by applying with a brush (brush-on technique).
 - Continue with Indirect method 9).
- Repairing temporary appliances
 - In case a repair is needed on a temporary appliance, remove the temporary appliance from the mouth and rough the fractured surfaces and adjacent areas. Undercuts should be provided. Apply mix material and allow the material to set. Details see above.

Repairs should not be done in the mouth.

B. REPAIR OF A FRACTURED DENTURE

- Take an impression of the edentulous area in a usual manner and fabricate a stone model.
- Roughen the fractured surfaces and adjacent areas of the denture using a bur and a fine diamond metal. The metal surface should be roughened with a diamond bur, sandblasted, cleaned with air and treated with METALPRIMER II.
- Coat the stone model with a separator such as ACROSEP or MULTI-SEP and place the fractured denture base on the stone model.
- Dispense the desired amount of UNIFAST III powder in a rubber cup and add liquid. Quickly mix with a plastic spatula for 10-15 seconds. Immediately pour the mixture into the fractured area of the denture base. Standard powder to liquid ratio is 1 g of powder to 0.5 mL of liquid (1 scale of powder) to 1 scale of liquid. Powder to liquid ratio is adjustable between 2 g : 1 mL and 2 g : 1.5 mL.
- Remove the model with pressure.
- Place the impression on the model with pressure.

- NOTE:** UNIFAST III starts to set 1 minute after mixing. Make sure to finish above steps within this period. If a longer working time is desired, store the powder and liquid in a refrigerator for more than 4 hours.
- Setting time is 3 minutes 10 seconds from start of mixing.
- After set, adjust, trim excess material and polish in the usual manner.

CLINICAL HINT – Brush on technique

As an alternative to mixing powder and liquid in a rubber cup, dispense the desired amount of powder and liquid directly on the fractured surface and adjacent areas. The liquid and pick up a small amount of powder. Deposit this mixture on the surface and progressively build up the temporary appliance, pattern or repair by applying the mixture on the fractured area of the denture base. Standard powder to liquid ratio is 1 g of powder to 0.5 mL of liquid (1 scale of powder) to 1 scale of liquid. Powder to liquid ratio is adjustable between 2 g : 1 mL and 2 g : 1.5 mL.

- NOTE:** UNIFAST III starts to set 1 minute after mixing. Make sure to finish above steps within this period. If a longer working time is desired, store the powder and liquid in a refrigerator for more than 4 hours.
- Setting time is 3 minutes 10 seconds from start of mixing.
- After set, adjust, trim excess material and polish in the usual manner.

IMPORTANT NOTE

When the uncured material is placed in the mouth, remove before the temperature starts to increase to avoid burn or necrosis of pulp.

- Do not use UNIFAST III to repair denture base materials or to reline pre-fabricated temporary crowns made out of polycarbonate, polysulfone, or polyethersulfone resin.
- Use of this product to repair metal surfaces is not recommended. If a metal surface of an appliance in water for at least 10 minutes prior to polishing in the patient's mouth.
- Small or thin portions of mixed material may not set properly. Make sure mixed material areas are fully cured before placing/plunging in the mouth.

SHADES

AI, A2, A3, A3.5, B2, E3, Incisal, Clear, No.3 (Pink), No.8 (Pink Veined)

PACKAGES

- UNIFAST III 2-Package
 - Kit contains 35 g of A3 (1), Powder 35 g A2 (1), Liquid 40 g (42 mL) (1), Plastic case (1), Liquid measure (1), Liquid measure (2), Disposable dispensing dish No.2 (5), Rubber cup (1), Plastic spatula (1), Brush No.10 (1)
- Refill Packages
 - Powder 35 g (in all shades)
 - Liquid 40 g (42 mL)
 - Powder 250 g (in all shades)
 - Powder 300 g (in all shades)
 - Liquid 100 g (104 mL)
 - Liquid 250 g (260 mL)

- Plastic Case
- Disposable dispensing dish No.2 (60)
- Brush No.10 (1)

*Only available in limited countries.

STORAGE

Recommended for optimal performance, store in a cool place (4-25°C, 39-77°F) away from direct sunlight.

CAUTION

- The liquid is flammable. Avoid sources of ignition. Do not keep or use near fire or high temperature.
- Avoid direct contact of powder, liquid or mixture on skin. In case of contact, wipe with a cotton pellet or sponge soaked in alcohol and flush with water.
- Avoid contact of powder, liquid or mixture on eyes. In case of contact, immediately flush with water and seek medical attention.
- The liquid is volatile. Use in a well-ventilated area. If mixed material is placed in the mouth before set, advise patients to breathe through the nose.
- Do not pour the powder of liquid remaining in the rubber cup directly into the bottle.
- Do not mix UNIFAST III with UNIFAST III or any other acrylic resin products.
- Do not use a brush or rubber cup which are used with other acrylic resin as this may delay or inhibit setting.
- Wear a dust protective mask when finishing and polishing the set material.
- Liquid may turn yellow with time, but this will not affect its function.
- The liquid should not be left in the rubber cup for prolonged periods of time as it changes the handling property.
- Personal protective equipment (PPE) such as gloves, face mask and safety eyewear should always be worn.
- In rare cases, the product may cause sensitivity in some people. If any such reactions do occur, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.georouse.com>


Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.georouse.com>

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for contamination. Disinfect immediately if damaged. Do NOT IMMERS. Thoroughly clean device to prevent drying and accumulation of contamination. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects: Reporting
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the links below:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vgilance@ge.dental
in this way you will contribute to improve the safety of this product.

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