Prior to use, carefully read the instructions for use. **EN**

GC Fuji™ I

RADIOPAQUE GLASS IONOMER LUTING CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Cementation of metal-based inlays, onlays, crowns, bridges and posts.
- Cementation of high strength (zirconia based) all ceramic crowns and bridges.

CONTRAINDICATIONS

- Pulp capping.
- Avoid use of this product in patients with known glass ionomer cement allergies.

DIRECTIONS FOR USE

Powder/liquid ratio (g/g)	1,8/1,0
Mixing time (min., sec.)	20"
Working Time (23°C / 73°F) (min., sec.)	2'00"
Finishing time (37°C / 99°F, after seating the restoration) (min., sec.)	4'30"

1. TOOTH PREPARATION

- Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
- Clean the prepared teeth with pumice and water.
- Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

2. RESTORATION PREPARATION

Make sure that the restoration is pretreated and handled according to the manufacturer’s instructions.

3. POWDER AND LIQUID DISPENSING

- The standard powder to liquid ratio is 1.8 g/1.0 g (1 level scoop of powder to 2 drops of liquid).
- For accurate dispensing of powder, lightly tap the bottle against the hand. DO NOT shake or invert.
- Hold the liquid bottle vertically and squeeze gently.
- Close bottles immediately after use.

4. MIXING

Dispense powder and liquid onto the pad. Using the plastic spatula, add all the powder to the liquid. Mix rapidly for 20 seconds.

Note:
When mixing larger amounts, divide the powder into two equal parts. Mix the first portion with all of liquid for 5 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15 seconds (total 20 seconds).

5. CEMENTATION

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes from start of mixing at 23°C (73°F). Higher temperatures shorten working time.
- Maintain moderate pressure.
- Start retaining excess cement when the excess cement feels rubbery.
- Finishing can be started 4 minutes 30 seconds after seating the restoration.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C) (39,2-77,0°F).

SHADE

Light yellow

PACKAGES

- 1-1 package: 35 g powder, 25 g (20 mL) liquid, powder scoop, mixing pads (No. 20).
- Bottle of 35 g powder with scoop.
- Bottle of 25 g (20 mL) liquid.

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with components of other glass ionomer cements.
- This product is not indicated for filling or core build-up.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

<http://www.gcasidental.com>
<http://www.gcaustralasia.com>
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

Last revised: 06/2019

MANUFACTURED BY
GC CORPORATION
78-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Rueschepark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL : +32 16 74 10 90

DISTRIBUTED BY
GC CORPORATION
78-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza, Singapore 529541
TEL: +65 6546 7598

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL :+61 2 9301 8200

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