Prior to use, carefully read the instructions for use.

EN

INLAY WAX SOFT INLAY WAX MEDIUM INLAY WAX HARD

DENTAL INLAY CASTING WAX

INLAY WAX SOFT, MEDIUM and HARD come in 3 different hardnesses and each has three different colors, so the most suitable wax can be selected to meet the needs of each case and technique. All these waxes exhibit ideal flow and plasticity, enabling wax build-up to be done with ease. They are hard to chip even in a thin formation. The wax patterns will not discolor working models even if they remain on the models for extended times. The wax containers are easy to handle with good function and design.

FEATURES

INLAY WAX SOFT

- Thanks to its excellent plasticity, INLAY WAX SOFT can easily be used for the build-up and pressure techniques.
- Because of its proper adaptability, INLAY WAX SOFT fits the die properly to prevent the wax from flaking at the margin due to its shrinkage.
- As INLAY WAX SOFT has the proper hardness for a soft type, it can easily be carved.

INLAY WAX MEDIUM

- The wax pattern will deform little at the time of solidification shrinkage, producing accurate restorations.
- INLAY WAX MEDIUM is designed specifically for the build-up technique. When it is softened, it exhibits excellent plasticity, so it can also be used for the pressure technique.
- The wax pattern is easy to remove from the model. There is no fear of deformation.
- As INLAY WAX MEDIUM has the proper hardness for a medium type, it can easily be carved.

INLAY WAX HARD

- INLAY WAX HARD can be carved soon after building up, providing for great efficiency in work.
- The wax is not sticky. Its chips will not stick to instruments and models.
- The wax pattern is easy to remove from the model. There is no fear of deformation.
- As INLAY WAX HARD has the proper hardness for a hard type, it can easily be carved.

PHYSICAL PROPERTIES

Test item	Solidification Temperature (°C)	Flow (%)			Thermal expansion (%) (25→37°C)
		30°C	40°C	45°C	
INLAY WAX SOFT	57.0	0.5	75.4	87.1	0.45
INLAY WAX MEDIUM	57.2	0.1	23.7	82.0	0.38
INLAY WAX HARD	60.6	0.1	1.1	56.9	0.40

APPLICATIONS

Wax up for inlays, crowns, bridges, etc.

DIRECTIONS FOR USE

- Apply a wax separating agent (such as GC SEP or Multi-Sep) to the surface of the model.
- If using the build-up technique, apply properly melted wax little by little using an instrument.
- If using the pressure technique, heat stick-shaped wax over a flame to adequately soften it, and apply the wax to the model.
- Perform carving and investing in the usual manner.
Note:
INLAY WAX SOFT, MEDIUM and HARD can be used with the GC Perfect Fit System which is available to take more accurate wax patterns.

NOTES

- INLAY WAX SOFT, MEDIUM and HARD are combustible. When using the material, take care not to get burnt.
- Do not leave the material at high temperatures (near a burning stove and under direct sunlight).
- Store in a cool, dark place away from fire. Do not keep large amounts in one single storage container.
- Fire extinguishing equipment should be installed in a place where the material is used or stored.
- These products should be used only by a dental professional.
- These products should be used only for the applications described in the instructions.

OTHER NOTES

- The inlay wax in gray contains a light shielding ingredient. It gives a different feeling of operation from waxes in red, violet or green of the same hardness.
- The same colors of three types may look different from one another.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<http://www.gceurope.com>
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:

MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

COLORS

- INLAY WAX SOFT: Gray, Red, Violet
- INLAY WAX MEDIUM: Gray, Red, Green
- INLAY WAX HARD: Gray, Red, Green

PACKAGES

- Tin of 40 g (SOFT, MEDIUM, HARD)
- Box of 60 sticks (SOFT, MEDIUM)

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GC ASIA DENTAL PTE. LTD.
GC AUSTRALASIA DENTAL PTY. LTD.
TAISEI DENTAL MFG. CO., LTD.