G-aenial™ Bond
ONE COMPONENT SELF-ETCHING LIGHT-CURED ADHESIVE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Bonding of light cured composites and acid modified composites (compomers) to tooth structures.
- Bonding of dual cured luting and core build up composites to tooth structure as long as these materials are light cured.

CONTRAINDICATIONS

- Pulp capping.
- Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer, or acetone.

NOT TO BE USED

- In combination with chemically cured composite resin.
- In combination with eugenol containing materials as eugenol may hinder G-aenial Bond from setting or bonding.
- In combination with desensitizers as desensitizers may hinder G-aenial Bond from setting or bonding properly.
- In combination with dual cure luting and core build up composites unless both G-aenial Bond and these materials are individually light cured.

DIRECTIONS FOR USE

1. CAVITY PREPARATION
Prepare tooth using standard techniques. Make sure that the temporary material is thoroughly removed from the surface. Use rubber dam for protection. Dry the prepared tooth surfaces by gently blowing with an air syringe.

2. TECHNIQUE SELECTION
Select out of the following 2 technique options
a) Self-etch technique – apply G-aenial Bond to bur cut enamel and dentin, without a separate etching step.
b) Selective etching of enamel – prior to applying G-aenial Bond to enamel and dentin, etch (bur cut) enamel with 35-40% phosphoric acid gel for 10 seconds, rinse for 5 seconds and gently dry.

3. APPLICATION

- Prior to dispensing, shake the bottle of G-aenial Bond thoroughly (Fig. 1). Dispense a few drops into a clean dispensing dish (Fig. 2). Replace bottle cap immediately after use (Fig. 3).
- IMMEDIATELY** apply to the prepared enamel and dentin surfaces using the dispensed applicator (Fig. 4).
- Leave undisturbed for 10 seconds after the end of application (Fig. 5). Then, dry thoroughly for 5 seconds with oil free air under MAXIMUM air pressure. Use vacuum suction to prevent splatter of the adhesive (Fig. 6). The final result should be a thin adhesive film with the appearance of frosted glass and which does not visibly move under further air pressure.

4. LIGHT CURING
Light cure using a visible light curing unit (Fig. 7).

IRADIATION TIME

Halogén / LED (700 mW/cm²)	: 10 segundos
High power LED (more than 1200 mW/cm²) :	5 segundos

In case where the light guide tip is more than 10mm away from the surface to be irradiated, light cure for the following time.

Halogén / LED (700 mW/cm²)	: 20 segundos
High power LED (more than 1200 mW/cm²) :	10 segundos

5A. APPLICATION OF LIGHT CURED COMPOSITES AND COMPOMERS
After having light cured the adhesive, follow manufacturer's instructions for use in order to place, contour and light cure the composite resin material of choice.

5B. PLACEMENT OF DUAL CURED COMPOSITES
After having light cured the adhesive, light cure the dual cured composite material separately. Self curing only will result in insufficient bonding strength.

6. FINISHING
Adjust restoration, finish and polish using standard techniques.

STORAGE
Store at room temperature (1-28°C) (33.8-82.4°F). If not in use for a prolonged period of time, store in refrigerator.

PACKAGE

- G-aenial Bond Bottle Kit: 5 mL liquid (1), DISPOSABLE DISPENSING DISH (20), DISPOSABLE APPLICATOR II (50)
- G-aenial Bond Bottle Refill: 5 mL liquid (1)
- G-aenial Bond 3-bottle Pack: 5 mL liquid (3)

CAUTION

- G-aenial Bond is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
- G-aenial Bond is volatile. Use in a well ventilated place. Advise patients to breath through nose.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
- If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks. To avoid contact, it is recommended to apply COCOA BUTTER to the area where rubber dam cannot cover.
- Avoid inhalation or digestion of material.
- Do not mix with other products.
- Dispose of all waste according to local regulations.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present FU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<http://www.gceurope.com>
or for The Americas
<http://www.gcamerica.com>

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects - Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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CE 2797

GC

MANUFACTURED BY
GC DENTAL PRODUCTS CORP.
2-285 Torinomiya-cho, Kasugagi, Aichi 486-8844, Japan
DISTRIBUTION IN NORTH AMERICA:
GC AMERICA INC.
3721 West 127th Street, Atlap, IL 60603 U.S.A.
DISTRIBUTION IN EUROPE:
GC CORPORATION
1-1-1, Higashi-Shinjuku, Shinjuku-ku, Tokyo 163-8585, Japan

GC SOUTH AMERICA
Rua Heitor Pires, 395, Santana, São Paulo, SP, BRASIL
CEP: 05022-051 - TEL: +55-11-2925-5065
FAX: +55-11-2925-5070
RESF. TEC. Mayara de Santa Ribeiro - C/RP/SP 105 982

GC ASIA DENTAL PTE. LTD.
5 Temasek Boulevard, #02-01, Temasek Place, Singapore 325541
TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1721 Boney Rd, Brisbane, QLD, NSW 4101, Australia
TEL: +61 2 9391 8200

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