




NEW METAL STRIPS

Colour-coded polishing strips for proximal surface



NEW METAL STRIPS

For preparation, polishing and finishing of proximal surface



Adequate flexibility with enough strength

A 18-8 stainless strip is carefully adopted to provide an adequate flexibility with enough strength to the push-and-pull movements.

Excellent durability

As fine particles of different sizes are electro-plated on 18-8 stainless strips, a longer durability is expected.

Repeated use with sterilization

Metal strips can be used repeatedly after sterilization either by medication or boiling. This may be recommended from the hygienic point.

Colour-coded

Colour code makes it easy to select the proper mesh.

Red mesh N° 200 - rough preparation

Blue mesh N° 300 - preparation

Green mesh N° 600 - polishing

Yellow mesh N° 1000 - finishing

Width - A(2.6mm), B(3.3mm), C(4.0mm)

Packages

12 pieces in a box

1. All same mesh - Box of 12 strips
(4 pieces each of 3 widths, A, B and C)
2. All same mesh - Box of 12 strips
(12 pieces - all one width)
3. Assortment - Box of 12 strips
(1 piece each width and each mesh)

GC CORPORATION

GC INTERNATIONAL CORP.

76-1Hassunuma-cho, Itabashi-ku, Tokyo 174, Japan Phone: 03-3558-5181 Fax: 03-3966-1470

For further information, please contact

GC BELGIUM, N.V. Research-Park, 13 Interleuvenlaan, B-3001 Leuven, Belgium Phone: 016-402700 Fax: 016-402214

GC GERMANY GmbH Robert-Bosch-Straße 17a 63477 Maintal-Dörmigheim, Germany Phone: 06181-947600 Fax: 06181-947629

GC FRANCE s. Avenue du Bouton d'Or, F-94380 Bonneuil s/Marne, France Phone: 01-4980-3791 Fax: 01-4980-3790

GC ITALY S.R.L. Via Montersato 62 (Sesto Uffiziario) I-20098 S. Giuliano Milanese (MI), Italy Phone: 02-96282068 Fax: 02-96282100

GC AMERICA INC. 3737 West 127th Street, Chicago, IL 60658 U.S.A. Phone: 708-597-0900 Fax: 708-371-5103

GC KOREA CO., LTD. #222, Changgeonno 3-Ga, Seodaejeon-Gu, Seoul, Korea Phone: 02-313-2272 Fax: 02-313-2275

GC TAIWAN CORP. 17 Lane 307, Chung-Shan Rd. Ser 2, Chung-Ho City, Taipei Hsien, Taiwan Phone: 2-246-3315 Fax: 2-245-5072

GC INTERNATIONAL CORP. SINGAPORE BRANCH

4th Floor, Jurdig Building, 421 Tagore Avenue, Singapore 3579 Phone: 65-4354700 Fax: 65-4366651

PRINTED IN JAPAN

160161/00