




COE® Disposable Trays

Plastic Impression Trays From GC-AMERICA

Spacer™ Trays,
Check-Bite™
Trays and
Coe Sideless
Triple Function
Trays



COE Brand



Spacer Trays

Full Arch - Perforated



- #1D Large Upper, 3" buccal width, 2-1/2" length, 7/8" anterior depth
- #4D Medium Upper, 2-7/8" buccal depth, 2-3/8" length, 3/4" anterior depth
- #7D Small Upper, 2-5/8" buccal width, 2" length, 11/16" anterior depth
- #20D Large Lower, 3" buccal width, 2-3/8" length, 3/4" anterior depth
- #21D Medium Lower, 3" buccal width, 2-1/4" length, 3/4" anterior depth
- #22D Small Lower, 2-15/16" buccal depth, 2" length, 13/16" anterior depth

Whether you need standard trays, triple function trays or bite registration trays, we have the economical, disposable and sturdy solution in the size you're looking for.

The unique zig-zag design of our **Spacer Trays** minimizes cusp contact with the tray and enables you to seat the tray to the correct depth for consistent, quality impressions. The anatomical design reduces pressure on soft tissue, and the strong, rigid construction provides excellent stability and support when removing impressions.

Disposable Impressio

The Broadest Selection





Full Arch - Solid

- #101D Large Upper, 3" buccal width, 2-1/2" length, 7/8" anterior depth
- #104D Medium Upper, 2-7/8" buccal depth, 2-3/8" length, 3/4" anterior depth
- #107D Small Upper, 2-5/8" buccal width, 2" length, 11/16" anterior depth
- #120D Large Lower, 3" buccal width, 2-3/8" length, 3/4" anterior depth
- #121D Medium Lower, 3" buccal width, 2-1/4" length, 3/4" anterior depth
- #122D Small Lower, 2-15/16" buccal depth, 2" length, 13/16" anterior depth

No matter what you need in disposable trays, plated metal trays or stainless steel trays, only GC America offers the COE brand, the widest selection of impression trays available.


**No one else
can make this statement.**

n Trays available anywhere.



DIMENSIONS

- A = Max. Buccal Width
- B = Max. Length
- C = Anterior Depth



Partial Arch / Quadrant Trays - Perforated



- # 30D Upper Left or Lower Right, 1-9/16" buccal width, 2-1/4" length, 3/4" anterior depth
- # 31D Upper Right or Lower Left, 1-9/16" buccal width, 2-1/4" length, 3/4" anterior depth
- # 32D Anterior, 2-9/16" buccal width, 1-3/8" length, 13/16" anterior depth

Partial Arch / Quadrant Trays - Solid

- #130D Upper Right or Lower Left, 1-9/16" buccal width, 2-1/4" length, 3/4" anterior depth
- #131D Upper Left or Lower Right, 1-9/16" buccal width, 2-1/4" length, 3/4" anterior depth
- #132D Anterior, 2-9/16" buccal width, 1-3/8" length, 13/16" anterior depth

GC®

Disposable Impression Trays




Check-Bite Triple Function

DIMENSIONS

A = Max. Buccal Width

B = Max. Length

C = Anterior Depth




Triple Function Trays

#73D Posterior Style, 2-3/8" buccal width, 2-3/16" length, 3/4" anterior depth

#75D Anterior Style, 2-5/8" buccal width, 1-1/2" length, 3/4" anterior depth

COE Sideless Triple Function



Sideless, Triple Function Trays

1st Hole 3" buccal width, 2-3/8" length, 1/8" anterior depth

2nd Hole 2-7/8" buccal width, 2-3/8" length, 1/8" anterior depth

3rd Hole 2-5/8" buccal width, 2-1/2" length, 1/8" anterior depth

4th Hole 2-1/2" buccal width, 2-1/2" length, 1/8" anterior depth

Triple Function Check-Bite Trays assure you of an accurate working impression and counter impression of the opposing arch, as well as a precise bite registration in one easy procedure. Check-Bite Trays cover more than half the arch area resulting in multiple bite registration points for more accurate articulation.

The sideless, **COE Triple Function Tray** is fully adjustable to 4 positions to accommodate various arch sizes and can be separated for taking a quadrant. Like our Check-Bite Trays, this versatile tray provides three impression functions in one easy step: working impression, counter arch impression and precise bite registration.



COE[®] Disposable Trays

Plastic Impression Trays

Description

Product Code

Spacer[™] Trays

Bag 12 Trays Shelf Pack 72 Trays Master Pack 432 Trays

Full Arch-Perforated

#1D Large Upper	250011	250016	250012
#4D Medium Upper	250041	250046	250412
#7D Small Upper	250071	250076	250712
#20D Large Lower	250201	250206	250212
#21D Medium Lower	250211	250216	252112
#22D Small Lower	250221	250226	252212
Assorted (Pkg. of each above)		250096	250112

Full Arch-Solid

#101D Large Upper	251011	251016	251012
#104D Medium Upper	251041	251046	251042
#107D Small Upper	251071	251076	251072
#120D Large Lower	251201	251206	251202
#121D Medium Lower	251211	251216	251212
#122D Small Lower	251221	251226	251222

Partial Arch / Quadrant Trays-Perforated

#30D Upper Left or Lower Right	250301	250306	250312
#31D Upper Right or Lower Left	250311	250316	253112
#32D Anterior, Perforated	250321	250326	253212

Partial Arch / Quadrant Trays-Solid

#130D Upper Right or Lower Left	251301	251306	251302
#131D Upper Left or Lower Right	251311	251316	251312
#132D Anterior	251321	251326	251322

Check-Bite[™]

Box

Triple Function Trays

#73D Posterior Style, box of 50	257350
#75D Anterior Style, box of 40	257540

COE Sideless, Triple Function

Box

Sideless, Triple Function Trays, box of 50 (25 Full Arch or 50 Quadrants)	150001
--	--------



75 Years of Advancing the Art and Science of Dentistry

GC AMERICA INC.

1-800-323-7063

1-800-GC FAX ME