

everStick®PERIO

Fibre reinforcement for periodontal splinting

DISTRIBUTED BY

GC CORPORATION

76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC EUROPE N.V.

Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium

TEL: +32 16 74 10 00

GC AMERICA INC.

3737 West 127th Street, Alsip, IL 60803 U.S.A., TEL: +1-708-597-0900

www.gcamerica.com

GC ASIA DENTAL PTE. LTD.

5 Tampines Central 1, #06-01 Tampines Plaza, Singapore 529541, TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.

1753 Botany Rd, Banksmeadow, NSW 2019, Australia, TEL: +61 2 9301 8200

GC SOUTH AMERICA

Rua Heliadora, 399, Santana - São Paulo, SP Brasil, TEL.: +55-11-2925-0965

CEP: 02022-051, CNPJ: 08.279.999/0001-61, RESP. TÈC: MAYARA DE SANTIS

RIBEIRO - CRO/SP 105.982

RESPONSIBLE MANUFACTURER IN CANADA

GC AMERICA INC.

3737 W. 127th Street, Alsip, IL 60803 U.S.A.



Fibre type: Silanated E-Glass fibre impregnated with Bis-GMA and PMMA

Form: Unidirectional fibre bundle

Diameter: ~ 1.1-1.3 mm



Keep away
from sunlight



Temperature
limit

951336 - NU6809
30000532

everStick®PERIO



- EN FIBRE REINFORCEMENT FOR PERIODONTAL SPLINTING
- DE GLASFASERVERSTÄRKUNG FÜR PARODONTALE SCHIENUNGEN
- FR FIBRE DE RENFORT POUR CONTENTION PARODONTALE
- IT RINFORZO IN FIBRA PER SPLINTAGGIO PARODONTALE
- ES FIBRA DE REFUERZO PARA FERULIZACIONES PERIODONTALES
- NL VEZELVERSTERKING VOOR PERIRADICULAIR SPALKEN
- DA FIBERFORSTÆRKET SPLINTING
- SV FIBERFÖRSTÄRKT SPLINTING
- PT FIBRA PARA REFORÇO DE CONTENÇÕES PERIODONTAIS
- EL ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΓΙΑ ΕΝΙΣΧΥΣΗ ΠΕΡΙΟΔΟΝΤΙΚΗΣ ΑΚΙΝΗΤΟΠΟΙΗΣΗΣ
- FI LASIKUITULUJITE KISKOTUKSIIN
- NO FIBERFORSTERKET SPLINTING
- ZH 复合树脂纤维用于牙周夹板治疗
- ID PENGUAT SERAT UNTUK SPLINTING PERIODONTAL

GC

Rx Only


