

everStick™ NET

Fibre reinforcement for labial splinting

DISTRIBUTED BY

GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium
TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A., TEL: +1-708-597-0900
www.gcamerica.com

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza, Singapore 529541, TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia, TEL: +61 2 9301 8200

GC SOUTH AMERICA
Rua Heliadora, 399, Santana - São Paulo, SP Brasil, TEL.: +55-11-2925-0965
CEP: 02022-051, CNPJ: 08.279.999/0001-61, RESP. TÊC: MAYARA DE SANTIS
RIBEIRO - CRO/SP 105.982

RESPONSIBLE MANUFACTURER IN CANADA

GC AMERICA INC.
3737 W. 127th Street, Alsip, IL 60803 U.S.A.

951334 - NU7406
30000530



Fibre type: Silanated E-Glass fibre impregnated with Bis-GMA and PMMA

Form: Bidirectional mesh fibre

Thickness: ~ 0,1 mm



Keep away
from sunlight



Temperature
limit

Caution: US Federal law restricts this device to sale by or on the order of a dentist." by **RxOnly**

everStick™ NET

- EN FIBRE REINFORCEMENT FOR LABIAL SPLINTING
- DE GLASFASERVERSTÄRKUNG FÜR LABIALE SCHIENUNGEN
- FR RENFORTS EN FIBRE DE VERRE POUR CONTENTION VESTIBULAIRE
- IT RINFORZO IN FIBRA PER SPLINTAGGIO LABIALE
- ES FIBRA DE REFUERZO PARA FERULIZACIONES LABIALES
- NL VEZELVERSTERKING VOOR LABIALE SPALKEN
- DA FIBERFORSTÆRKNING AF TRAUMER
- SV FIBERFÖRSTÄRKNING FÖR LABIAL SPLINTING
- PT FIBRA DE REFORÇO PARA CONTENÇÃO VESTIBULAR
- EL ΕΝΙΣΧΥΣΗ ΓΛΩΣΣΙΚΗΣ ΑΚΙΝΗΤΟΠΟΙΗΣΗΣ ΜΕ ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ
- FI LASIKUITULUJITE LABIAALISEEN KISKOTUKSEEN
- ZH 复合树脂纤维 用于唇牙周夹板治疗
- ID SERAT PENGUAT UNTUK SPLINTING LABIAL

GC



